

GLOBAL MEDICAL PRODUCTS

IN CANADA: 1-800-387-6095

AU QUÉBEC: 1-800-655-9445

info@globalmedical.ca



PEDIATRIC SCALES



**SCALE-TRONIX
4802 Pediatric Scale
4302 Diaper /
Specimen
Scale**

4802

For Weighing Neonates through Pediatrics

- Automatic Zero
- Weight recall.
- Reweigh.
- Absolute weight, no averaging.
- Accurate to 1 gm/5 gm.
- Optional Diaper weighing & Weight printer
- Disposable liners for infection control.
- 45 lb/20 kg capacity.
- Bright, easy-to-read digital display.

The 4802 is designed for daily and serial weighing of neonates, infants and toddlers. It's a heavy-duty, stainless steel, portable unit with a built-in digital readout. This scale is accurate enough to weigh neonates and has a large enough capacity for peds. The 4802 automatically returns to zero, allowing you to weigh a continuous stream of patients quickly and efficiently. A "Recall" feature retains the last weight reading. If the patient is removed from the scale and you forget the weight, simply press the "Recall" button. The previous weight is displayed, and a "Prior Weight" message blinks.

Since accurate weighing of active babies is difficult, many scales take an average between a high and low reading or lock-in a reading after an arbitrary length of time. Scale-Tronix systems display the actual weight so you will never medicate or treat an infant based on an inaccurate reading. In addition, the press of a button allows reweighing while the baby is still on the scale. An audible signal indicates when the patient's weight is stored in memory. A printer is available to provide a paper-tape record of patient weights.

The 4802's readings are so accurate and repeatable that breast-feeding intake may be determined by weight. The 1 gram weighing option enables you to weigh diapers automatically to 1 gram, ideal for determining urine output.

The 4802 is line-cord powered for stationary use or powered by its self-contained battery for portable use. Weight is displayed with the touch of a button in either pounds or kilos. Oversized and 4-sided cradles are also available. A built-in dispenser holds extra-large disposable liners used for infection control.

The unique Scale-Tronix stainless steel weighing cart makes the 4802 a self-contained, mobile weighing station. The cart also features a dispenser for disposable scale liners and two full-size shelves. Computer output capability can be added.

MORE THAN 700 IN USE IN CANADA



Specifications: 4802 Pediatric Scale

ACCURACY: 2/10 OZ/5 GM. DIAPER/PREEMIE WEIGHING 1 GM OPTION).

WEIGHING RANGE: 45 LB/20 KG.

PLATFORM WITH CARRYING HANDLES:

OVERALL, 27 1/2" X 16" X 7" STAINLESS STEEL.

REMOVABLE CRADLE: ACRYLIC WITH LENGTH MEASURE IN IN/CM.

25" X 14", STANDARD. ALSO: OVERSIZED 32" X 14",

FOUR SIDED 25" x 14" (Available as options)

READOUT: DIGITAL LED DISPLAY IN POUNDS AND KILOGRAMS.

POWER SOURCE: 110 VAC LINE CORD AND SELF-CONTAINED BATTERY WITH CHARGER (OPERATES IN EITHER MODE).

AUTOMATIC ZERO: RETURNS TO ZERO, READY FOR NEXT PATIENT.

REWEIGH: RECOMPUTES WEIGHT WHILE PATIENT IS ON SCALE.

WEIGHT RECALL: RECALLS LAST WEIGHT.

LOW BATTERY INDICATOR.

CHARGE INDICATOR: INDICATES BATTERY IS CHARGING.

WEIGHING CART ON WHEELS: STAINLESS STEEL, WITH DISPENSER FOR CRADLE LINERS AND TWO STORAGE SHELVES. 33"h X 21 1/2"d X 32 1/2"w. PRINTER: ACCESSORY.

COMPUTER CAPABILITY: RS 232 OUTPUT (OPTIONAL).
UL/CSA LISTED.

MODEL 4802 OPTIONS / ACCESSORIES:

48045, 1 GM ACCURACY.

48800-1 MOBILE WEIGHING CART WITH
DISPENSER FOR SCALE LINERS & STORAGE SHELVES

48220 FOUR-SIDED CRADLE*.

48222 OVERSIZED CRADLE*.

* WHEN ORDERED IN PLACE OF STANDARD CRADLE



50% OF CUSTOMERS PLACE THEIR SCALE ON THE MOBILE CART.

4302

Organ/Specimen Diaper Scale



- Accurate to 1 gm.
- Push-Button zero.
- Battery power.
- Bright, easy-to-read digital display.
- Overload protection.

The 4302 was created in consultation with surgeons and nurses to ensure suitability to its crucial task in O.R. pathology and general hospital use.

Made of durable, easy-to-clean stainless steel, this compact unit has multiple handles for ease of use and transport. The 8 1/2" x 8 1/2" platform accommodates all typical bowls or basins.

Specifications: 4302 Organ/Specimen Diaper Scale

ACCURACY: 1 GM.
WEIGHING RANGE: 0-3,000 GM.
READOUT: DIGITAL 1/2" (RED) LED
IN GRAMS, POUNDS OR OUNCES.
AUTOMATIC SHUTOFF.
PLATFORM SIZE: 8 1/2" X 8 1/2".
POWER SOURCE: 6 "C" SIZE
DISPOSABLE ALKALINE BATTERIES.
EXTERNAL ADAPTER IS AVAILABLE
FOR LINE-POWER OPERATION.

The low-profile scale's proprietary weighing mechanism provides rugged, yet sensitive operation. One gram (.0001 lb.) resolution assures precision weighing of even small items. Capacity is 3,000 gm (6.6 lb.) with overload protection built in. Weight is displayed on an easy-to-read digital readout. With the push of a button, the scale zeros out the container weight providing net weight of the contents.

Six self-contained, size "C" alkaline batteries provide reliable operation

BATTERY LIFE: IN EXCESS OF 10,000
WEIGHINGS.
LOW BATTERY INDICATOR.
CONSTRUCTION: STAINLESS STEEL.
HANDLES: 3 HANDLES PROVIDED.
DIMENSIONS: 8" w X 11" d X 3" h.
WEIGHT: 8 LBS.



for more than a year. A time-out feature automatically turns the scale off after use, extending battery life. Users are alerted by a low-battery indicator well before batteries fade. For continuous operation, an AC power adapter is optional.

Its sensitivity and robust construction make the 4302 ideal for a number of other applications: weighing laboratory animals, sanitary pads for mothers' blood loss, breast reduction, O.R. sponges, etc.

MODEL 4302
OPTION:
EXPANDED WEIGHING CAPACITY
TO 5,000 GRAMS.
341001 POWER ADAPTER,
070070 EXTRA MANUAL.